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Tıbbi Cihaz

Eng3 Şirketi

Türk

M004-TUR-rev19

**eng3**

Life Science Technology

eng3

# NanoVi® Eco / NanoVi® Pro / NanoVi® Exo

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eng3

## İçindekiler

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<b>Guidance and manufacturer's declaration - electromagnetic emissions</b>		
The NanoVi is intended for use in the electromagnetic environment specified below. The customer or the user of the NanoVi should assure that it is used in such an environment.		
<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment – guidance</b>
RF emissions CISPR 11	Group 1	The NanoVi uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The NanoVi is suitable for use in all establishments, including domestic and medical establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

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<b>Guidance and manufacturer's declaration - electromagnetic immunity</b>			
The NanoVi is intended for use in the electromagnetic environment specified below. The customer or the user of the NanoVi should assure that it is used in such an environment.			
<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment - guidance</b>
Electrostatic discharge (ESD)  IEC 61000-4-2	± 8 kV contact ± 2 kV air ± 4 kV air ± 8 kV air ± 15 kV air	± 8 kV contact ± 2 kV air ± 4 kV air ± 8 kV air ± 15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %
Electrical fast transient/burst  IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for SIP/SOP Repetition frequency 100 kHz	± 2 kV for power supply lines Repetition frequency 100 kHz	Mains power quality should be that of a typical commercial or hospital environment.
Surge  IEC 61000-4-5	± 0,5, 1,0 kV line(s) to line(s) ± 0,5, 1,0, 2 kV line(s) to earth	± 0,5, 1,0 kV line(s) to line(s) ± 0,5, 1,0, 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.

<p>Voltage dips, short interruptions and voltage variations on power supply input lines</p> <p>(Blackouts, brownouts, and fluctuations of the power supply according to IEC )</p> <p>IEC 61000-4-11</p>	<p>0 % UT; 0.5 cycle</p> <p>At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°</p> <p>0 % UT; 1 cycle and 70 % UT; 25/30 cycles</p> <p>Single phase: at 0°</p> <p>0 % UT; 250/300 cycles</p>	<p>0 % UT; 0.5 cycle</p> <p>At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°</p> <p>0 % UT; 1 cycle and 70 % UT; 25/30 cycles</p> <p>Single phase: at 0°</p> <p>0 % UT; 250/300 cycles</p>	<p>Mains power quality should be that of a typical commercial or hospital environment. If the user of the NanoVi requires continued operation during power mains interruptions, it is recommended that the NanoVi be powered from an uninterruptible supply or a battery.</p>
<p>Power frequency (50/60Hz) magnetic field</p> <p>IEC 61000-4-8</p>	<p>30 A/m</p>	<p>30 A/m</p>	<p>Power frequency magnetic fields should be at levels characteristic of a typical commercial or hospital environment.</p>
<p>NOTE U<sub>T</sub> is the a.c. mains voltage prior to application of the test level.</p>			

**Guidance and manufacturer's declaration - electromagnetic immunity**

The NanoVi is intended for use in the electromagnetic environment specified below. The customer or the user of the NanoVi should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
<p>Conducted FR</p> <p>IEC 61000-4-6</p>	<p>3 Vrms</p> <p>150 kHz to 80 MHz</p> <p>6 Vrms in ISM bands between 0.15 MHz and 80 MHz</p> <p>80 % AM at 1 kHz</p>	<p>3 V</p> <p>150 kHz to 80 MHz</p> <p>6 Vrms in ISM bands between 0.15 MHz and 80 MHz</p> <p>80 % AM at 1 kHz</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the NanoVi, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are normally.</p> <p><b>Recommended separation distance</b></p> $dd = \frac{3.5}{3} \sqrt{PP}$ <p><math>dd = \frac{3.5}{10} \sqrt{PP}</math>      80 MHz to 800 MHz</p> <p><math>dd = \frac{7}{10} \sqrt{PP}</math>      800 MHz to 2,5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p>
<p>Radiated RF</p> <p>Transient RF</p> <p>IEC 61000-4-3</p>	<p>3 V/m</p> <p>80 MHz to 2,7 GHz</p> <p>80 % AM by 1 kHz</p>	<p>10 V/m</p> <p>80 MHz to 2,7 GHz</p> <p>80 % AM by 1 kHz</p>	

			<p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency band.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol.</p> 
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structure, objects, and people.</p> <p><sup>a</sup> Field strengths from fixed transmitters, such as base stations for radios (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the NanoVi is used exceeds the applicable RF compliance level above, the NanoVi should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the NanoVi.</p> <p><sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

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Recommended separation distances between portable and mobile RF communications equipment and the NanoVi			
<p>The NanoVi is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the NanoVi can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the NanoVi as recommended below, according to the maximum output power of the communications equipment.</p>			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $dd = \frac{3,5}{3} \sqrt{PP}$	80 MHz to 800 MHz $dd = \frac{3,5}{10} \sqrt{PP}$	800 MHz to 2,5 GHz $dd = \frac{7}{10} \sqrt{PP}$
0,01	0,12	0,035	0,07
0,1	0,37	0,11	0,22
1	1,2	0,35	0,7
10	3,8	1,1	2,2
100	12	3,5	7
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.</p> <p>NOTE 2 The guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

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	NanoVi® Eco	NanoVi® Pro	NanoVi® Exo
Article Number:	4650-00	4800-00	4900-00
Output performance: NanoVi® Eco has half the output: 50%  NanoVi® Pro establishes baseline output: 100%  NanoVi® Exo has 2x the output: 200%			
Recommended frequency of sessions:	Optimal: Daily Ideal: 3 x per week Minimum: 1 x per week	Optimal: Daily Ideal: 3 x per week Minimum: 1 x per week	Optimal: Daily Ideal: 3 x per week Minimum: 1 x per week
Suggested session time (More is better)	<b>60 minutes on the Eco</b> =	<b>30 minutes on the Pro</b> =	<b>15 minutes on the Exo</b>
Minimum lifespan:	10,000 hours = 600,000 minutes	10,000 hours = 600,000 minutes	10,000 hours = 600,000 minutes
Minimum # of sessions per lifespan:	<b>10,000 sessions of 60 minutes</b>	<b>20,000 sessions of 30 minutes</b>	<b>40,000 sessions of 15 minutes</b>
Smartcard system: - <b>Owner Card</b> - <b>User Card</b> works on the device it was programmed for	No	Yes: - Owner Card for selecting the Operation Mode - User Card programmable for up to 9,999 minutes	Yes: - Owner Card for selecting the Operation Mode - User Card programmable for up to 9,999 minutes
Standard mode - Timer Mode:	<b>Yes</b> , enter application time in minutes on the touch screen.	<b>Yes</b> , enter application time in minutes on the touch screen.	<b>Yes</b> , enter application time in minutes on the touch screen.
Mode with Smartcard - Timer Mode:	<b>No</b>	<b>Yes</b> , used application time gets deducted from the User Card.	<b>Yes</b> , used application time gets deducted from the User Card.
- Session Mode:	<b>No</b>	<b>Yes</b> , session minutes deducted from User Card.	<b>Yes</b> , session minutes deducted from User Card.
- Infinity Mode:	<b>No</b>	<b>Yes</b>	<b>Yes</b>
Option for grouping multiple devices: (With Smartcard System only)	<b>No</b>	<b>Yes</b> , User Card can be used for any device in the group.	<b>Yes</b> , User Card can be used for any device in the group.
Lamp for illumination of water: Illumination effects:	16 Colors & off Solid, flash, strobe, fade, smooth	16 Colors & off Solid, flash, strobe, fade, smooth	16 Colors & off Solid, flash, strobe, fade, smooth
Excitation units / - elements / max power:	1 / 6 / 6 x 700 pW	2 / 12 / 12 x 700 pW	4 / 24 / 24 x 700 pW
Spectral emission / max power:	1100–1300nm and 1500-1700nm	1100–1300nm and 1500-1700nm	1100–1300nm and 1500-1700nm
Dimensions (w x l x h): Weight:	12" x 11" x 9" (31 x 28 x 23 cm) 8.5 lb. (3.9 kg)	12" x 11" x 9" (31 x 28 x 23 cm) 9.0 lb. (4.1 kg)	12" x 11" x 9" (31 x 28 x 23 cm) 9.5 lb. (4.3 kg)
Silver antimicrobial tubing: Amount of distilled water for operation:	Yes, after humidification unit 8.5 fl. oz. (250 ml)	Yes, after humidification unit 8.5 fl. oz. (250 ml)	Yes, after humidification unit 8.5 fl. oz. (250 ml)
Display: Volume and brightness adjustment:	Touch Screen LCD, color 1 factory setting	Touch Screen LCD, color 8 individual settings	Touch Screen LCD, color 8 individual settings
Pumps: Air intake and output:	2 0.141 CFM (4.0 liters per minute)	2 0.141 CFM (4.0 liters per minute)	2 0.141 CFM (4.0 liters per minute)
Current Levels (A) Off	0.04	0.04	0.04
Current Levels (A) On-standby	0.10	0.10	0.10
Current Levels (A) On-running	0.30	0.32	0.36
Components automatically tested for function:	Pumps, excitation elements, lamp	Pumps, excitation elements, lamp	Pumps, excitation elements, lamp
Service check:	Not Required, self-testing	Not Required, self-testing	Not Required, self-testing
Warranty:	2 Years	2 Years	2 Years

#### Accessories

4300-00	Glass Container	Glass
4370-00	Holder with Glass Container	Aluminum
4340-00	Power Supply Input / Output	100 - 240V AC / 12V DC, 4A
4350-00	Power Cord 2-Prong, 6 ½ ft. (2 m)	US, CA, Mexico, Japan
-10	Power Cord 2-Prong, 6 ½ ft. (2 m)	EU, Asia, Israel, South America
-20	Power Cord 2-Prong, 6 ½ ft. (2 m)	UK, Malaysia, Singapore
-30	Power Cord 2-Prong, 6 ½ ft. (2 m)	Australia, New Zealand
-60	Power Cord 2-Prong, 6 ½ ft. (2 m)	China
4530-00	User Manual	Available in different languages
4540-00	Concise User Manual	Available in different languages
<b>Consumables</b>		
4310-00	Diffuser for Humidifier	Replace once a year
4400-00	Paper tube	White paper tubes

#### Optional Accessories

4200-00	Transportation Case	Black with foam inserts
<b>Consumables</b>		
4360-00	Nasal cannulas	1 ft. (0.3 m), Latex-free material
4510-00	User Smartcards for NanoVi® Pro and NanoVi® Exo	Programmable for up to 9,999 minutes for use in "Card Mode"
4510-00	Owner Smartcards for NanoVi® Pro and NanoVi® Exo	Required for setting "Card Mode" and for programming minutes on User Smartcards

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U.S. Department of Health and Human Services Welcome, Hans-Joachim - ENG | FURLS HOME

**FDA FURLS** **DRLM**  
 Device Registration & Listing Module

DRLM Home > View Your Registrations and Listings

**Annual Registration** **View Your Registered Facilities**  
 Annual Registration  
 Owner/Operator: 9058276

**Facility Registration**  
 Register a New Medical Device Facility  
 Change Registration Information for a Facility  
 Cancel, Deactivate, or Reactivate a Facility Registration  
 View Your Registration and Listing Information

Show 10 per page [Clear Sort and Filter](#)  
Filter:

Registration Number	Registration Status	Registration Status Reason	Facility Name/ Trade Name(s)	Address	Expiration Date	Action
3004152208	Active	Registration changed from inactive to active	ENG3 CORPORATION	2234 EASTLAKE AVE E. SEATTLE, WASHINGTON 98102 UNITED STATES	2025-12-31	<a href="#">View</a>

DRLM Home > View Your Registrations and Listings

**Annual Registration** **View Your Device Listings**  
 Annual Registration  
 Owner/Operator: 9058276

**Facility Registration**  
 Register a New Medical Device Facility  
 Change Registration Information for a Facility  
 Cancel, Deactivate, or Reactivate a Facility Registration  
 View Your Registration and Listing Information  
 Change the Official Correspondent

Show 10 per page [Clear Sort and Filter](#)  
Filter:

Listing Number	Listing Status	Premarket Submission Number	Product Code(s)	Type of Combination Product	Device Name	Registration Number/FEI [Activities]	Action
D097353	Active		KFZ		HUMIDIFIER, NON-DIRECT PATIENT INTERFACE (HOME-USE)	Registration Number: 3004152208/3004152208 [Manufacturer]	<a href="#">View</a>

DRLM Home > View Your Registrations and Listings

**Annual Registration** **View Selected Listing Details**  
 Annual Registration  
 Listing Number: D097353  
 Listing Status: Active  
 Premarket Submission Number:

**Facility Registration**  
 Register a New Medical Device Facility  
 Change Registration Information for a Facility  
 Cancel, Deactivate, or Reactivate a Facility Registration  
 View Your Registration and Listing Information  
 Change the Official Correspondent for a Facility

**Facility Ownership**  
 Transfer Ownership of a Facility (Reopt Purchase)

Product Code	Product Name
KFZ	HUMIDIFIER, NON-DIRECT PATIENT INTERFACE (HOME-USE)

Registration #	Registration Status	Registration Status Reason	Activities
3004152208	Active	Registration changed from inactive to active	Manufacturer

[View Proprietary Names](#)

DRLM Home > View Your Registrations and Listings

**View Proprietary Names**  
 Listing Number: D097353

Proprietary Name	Confidential
NanoVi	N
NanoVi Exo	N
NanoVi Pro	N
NanoVi Eco	N

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2797

Salter Labs  
30 Spur Drive  
El Paso, Texas 79906  
USA

**DECLARATION OF CONFORMITY**  
DOC-20017a Oxygen Cannulas

SRN: US-MF-00007934  
Basic UDI: 0060741F1005\_Cannulas27

5Product Name:  
Intended Purpose

Model Number or Designator:

Control Designator:

Device Classification:

Conformity Route:

EMDN Nomenclature code (EMDN):

Global Medical Device Nomenclature Code (GMDN):

Universal Medical Device Nomenclature System (UMDNS):

Product Options/Accessories:

EC Certificate – Full Quality Assurance System

Notified Body:

Authorized EU Representative:

Revision: 1

Authorized Signature and Fund  
Printed Name: Rob Yamashita,  
Vice President Regulatory

Date of Issue: 12-Sept-2022  
Place of Issue: Regulatory Office, SunMed  
Grand Rapids, Michigan, USA

This declaration of conformity is issued under the sole responsibility of SunMed. We hereby declare that the medical device(s) specified above meet the provision of the Regulation (EU) MDR 2017/745 for medical devices.  
All technical documentation is retained at the premises of the manufacturer/technical documentation location.

Doc #: FRM-10002 Rev 3 Page 1 of 12

2797

Salter Labs  
30 Spur Drive  
El Paso, Texas 79906  
USA

**DECLARATION OF CONFORMITY**  
DOC-20017a Oxygen Cannulas

SRN: US-MF-00007934  
Basic UDI: 0060741F1005\_Cannulas27

UDI # GTIN	Model Number / Case Qty	Description
00607411919145	13320	Nasal Cannula (Premature) Salter Style with 7' (2.1 m) supply tube.
30607411919146	13320-50	
00607411919152	13321	Nasal Cannula (Neonate) Salter Style with 7' (2.1 m) supply tube.
30607411919153	13321-50	
00607411919169	13322	Nasal Cannula (Infant) Salter Style with 7' (2.1 m) supply tube.
30607411919160	13322-50	
00607411919176	13323	Nasal Cannula (Intermediate Infant) Salter Style with 7' (2.1 m) supply tube.
30607411919177	13323-50	
00607411919183	13324	Nasal Cannula (Pediatric) Salter Style with 7' (2.1 m) supply tube.
30607411919184	13324-50	
00607411100000	1600-1	Nasal Cannula (adult) Salter Style® with 1' (0.3 m) supply tube.
10607411100007	1600-1-50	
00607411100017	1600-10	Nasal Cannula (adult) Salter Style® with 10' (3.0 m) supply tube.
10607411100014	1600-10-50	
00607411100031	1600-12	Nasal Cannula (adult) Salter Style® with 12' (3.65 m) supply tube.
10607411100038	1600-12-50	
00607411000048	1600-13	Nasal Cannula (adult) Salter Style® with 13' (4.0 m) supply tube.
10607411100499	1600-13-50	
00607411100055	1600-14	Nasal Cannula (adult) Salter Style® with 14' (4.3 m) supply tube.
10607411100052	1600-14-50	
00607411100062	1600-15	Nasal Cannula (adult) Salter Style® with 15' (4.6 m) supply tube.
10607411100069	1600-15-50	
00607411100079	1600-16	Nasal Cannula (adult) Salter Style® with 1' (0.3 m) supply tube.
10607411100076	1600-16-50	
00607411100109	1600-2	Nasal Cannula (adult) Salter Style® with 2' (0.6 m) supply tube.
10607411100106	1600-2-50	
00607411100116	1600-20	Nasal Cannula (adult) Salter Style® with 20' (6.1 m) supply tube.
10607411100113	1600-20-50	
00607411100123	1600-21	Nasal Cannula (adult) Salter Style® with 21' (6.4 m) supply tube.
10607411100122	1600-21-50	
00607411100147	1600-24	Nasal Cannula (adult) Salter Style® with 24' (7.3 m) supply tube.
10607411100144	1600-24-50	
00607411100154	1600-25	Nasal Cannula (adult) Salter Style® with 25' (7.6 m) supply tube.
10607411100151	1600-25-50	
00607411100161	1600-3	Nasal Cannula (adult) Salter Style® with 3' (0.9 m) supply tube.
10607411100168	1600-3-50	
00607411100178	1600-30	Nasal Cannula (adult) Salter Style® with 30' (9.1 m) supply tube.
10607411100175	1600-30-20	
00607411100185	1600-35	Nasal Cannula (adult) Salter Style® with 35' (10.7 m) supply tube.
10607411100182	1600-35-20	

Doc #: FRM-10002 Rev 3 Page 3 of 12

17.2.4 NanoVi® Kanül Kaydı: FDA

U.S. FOOD & DRUG ADMINISTRATION

Home Food Drugs Medical Devices Radiation-Emitting Products Vaccines, Blood & Biologics Animal & Veterinary Cosmetics Tobacco Products

### Establishment Registration & Device Listing

FDA Home Medical Devices Databases

New Search Back To Search Results

**Establishment:**  
 Salter Labs  
 30 Spur Drive  
 El Paso, TX 79906  
**Registration Number:** 3004748541  
**FEI Number\*:** 3004748541  
**Status:** Active  
**Initial Distributor/Importer:** Yes  
 \*Note Firm May Have Additional Establishment Types.  
 Please Review Listings For Further Information.  
**Date Of Registration Status:** 2025

**Owner/Operator:**  
 SunMed  
 2710 Northridge Drive NW, Suite A  
 Grand Rapids, MI US 49544  
**Owner/Operator Number:** 9038422

**Official Correspondent:**  
 Rob Yamashita  
 2710 Northridge Drive NW, Suite A  
 Grand Rapids, MI 49544  
**Phone:** 01-616-2598373

\* Firm Establishment Identifier (FEI) should be used for identification of entities within the imports message set

Page Last Updated: 01/06/2025

U.S. FOOD & DRUG ADMINISTRATION

Home Food Drugs Medical Devices Radiation-Emitting Products Vaccines, Blood & Biologics Animal & Veterinary Cosmetics Tobacco Products

### Establishment Registration & Device Listing

FDA Home Medical Devices Databases

New Search Back To Search Results

**Proprietary Name:** Salter Labs nasal cannulas, 16Soft series, 16 series, high flow series, OTC series, gator and ear wraps accessory

**Classification Name:** CANNULA, NASAL, OXYGEN

**Product Code:** CAT

**Device Class:** 1

**Regulation Number:** 868.5340

**Medical Specialty:** Anesthesiology

**Registered Establishment Name:** Salter Labs

**Registered Establishment Number:** 3004748541

**Owner/Operator:** SunMed

**Owner/Operator Number:** 9038422

**Establishment Operations:** Specification Developer

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CH REP ALBO-Healthcare GmbH  
Alte Steinhäuserstrasse 19  
6330 Cham, Switzerland  
+41.43.818.75.02

UK REP UK RepMed Ltd  
Coxbridge Business Park, Unit D,  
Crosdail Place, Alton Rd,  
Farnham GU10 5EH, UK  
+44.1252.912.933

Kalms Consulting GmbH  
Rheinstr. 45-46  
12161 Berlin, Germany  
+49.30.405.045.320

↑  
Yetkili  
Temsilci

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Eng3 Corporation  
2234 Eastlake Ave E  
Seattle, WA, USA  
[www.eng3corp.com](http://www.eng3corp.com)

[www.eng3corp.com/service](http://www.eng3corp.com/service)

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
**Nasal Cannula** Home Oxygen Instructions for Use  
Please read and follow the Instructions for Use prior to using your nasal cannula for your home oxygen therapy.

The Nasal Cannula is used to deliver supplemental oxygen to patients who have a prescription for home oxygen therapy. Prior to home use, you and/or your caregiver should receive instructions from a trained healthcare professional on how to safely use your nasal cannula while on oxygen.

Nasal Cannulas are disposable and for single-patient use.

Nasal Cannulas are intended for use in the home, outpatient, extended care, transport and hospital environments.

Nasal cannulas are available in sizes from infants to adults. The nasal cannulas are available in different styles with various lengths of supply tubing. Some styles may have liter flow limitations, which will be stated on the product label.

Physical problems (e.g., chest pains, cannot breath), call 911.

#### Instructions for Use (continued)

##### Safety Precautions

- Oxygen is a nonflammable gas, but does support combustion. Follow your homecare provider's instructions for the care and safe operation of your oxygen delivery system (e.g., oxygen cylinder, oxygen concentrator, liquid oxygen).
- Do not smoke or allow anyone to smoke around you.** This includes, but is limited to, cigarettes, pipes, cigars, and electronic cigarettes (vapers). 
- Keep oxygen equipment at least 6 feet away from flames** or any heat source, for example, fireplaces, stoves, barbeque grills, and space heaters. 
- Example, don't apply Vaseline around or in your nares. 
- Do not use flammable products such as aerosol sprays or cleaning products while wearing your nasal cannula or around your oxygen source. 
- Avoid using electrical equipment that may cause a spark, for example, electric razor, blow-dryer or curling iron. 
- Use oxygen as prescribed by your doctor.
- The total length of your nasal cannula and oxygen supply tubing should not exceed 57 feet to ensure there is enough pressure to deliver prescribed oxygen flow rate.
- Do not kink, bend or tie your oxygen tubing.
- Do not place anything on your tubing that may obstruct flow.
- Keep excess tubing loosely coiled and out of the way to prevent tripping on oxygen tubing.
- Do not place your oxygen tubing or nasal cannula under blankets, bedsheets, rugs, etc.
- Use caution to prevent your oxygen tubing from becoming entangled in your furniture.
- Keep an extra nasal cannula and other oxygen supplies available for use.
- Do not let children or pets play with your nasal cannula and oxygen equipment.
- Recommend use of swivel adapter for nasal cannula and supply tubing longer than 14 feet.
- If using humidification, add a water trap to collect excess moisture in the supply tubing.



(Continued on the next page)

#### Instructions for Use (continued)

Check for gas flow from the nasal prongs. low.

5. a—Wrap the headset loop up and over both ears.

6. a—Squeeze the sides of the bolo and glide the bolo up under your chin.  
b—Leave enough space to fit at least two fingers between the bolo and chin.

(Continued on the next page)

SLML-130 Rev B, Aug 2016

#### Instructions for Use (continued)

##### Troubleshooting Tips

Problem	Possible Cause	Corrective Action
No oxygen flow from nasal prongs	1. Cannot feel the airflow in your nostrils. 2. Flow control valve is not turned on. 3. Oxygen system is not functioning properly or oxygen container is empty. 4. The nasal cannula is disconnected from oxygen device or supply tubing. 5. Nasal cannula or oxygen tubing kinked or blocked.	1. Check air flow by placing prongs next to hand or place nasal prongs into a small container of clean water. Bubbles will appear if there is oxygen flow. 2. Set flow control to prescribed setting. 3. Switch to backup oxygen source and contact your homecare provider. 4. Reconnect oxygen tubing. Ensure all tubing connections are tight and secure. 5. Inspect cannula and oxygen tubing for kinks or damage. Ensure nothing is placed on top of the tubing.
Water in nasal cannula or oxygen supply tubing	1. Humidifier bottle overfilled, or bottle has tipped over. 2. Water trap is full 3. High humidity environment, or sudden drop in temperature.	1. Pour out the excess water. Ensure that the humidifier bottle is upright. 2. Empty water trap. 3. Consider adding a water trap to your oxygen supply tubing.
Nasal dryness or irritation	1. Gas flow is dry. 2. No humidifier is being used.	1. Use normal saline spray or water soluble ointment, (i.e. AYR Saline Nasal Gel) to moisten the inside your nostrils. If condition worsens, contact your doctor. 2. Contact your doctor or homecare provider to request humidification.
Soreness or irritation around ears	1. Headset tubing too tight. 2. Tubing pressing against skin.	1. Loosen headset tubing. 2. Place a cotton padding or cushion (i.e., EZ-Wrap) under headset tubing.
Skin rash and/or sores caused by nasal cannula	1. Sensitivity or reaction to nasal cannula material. 2. Nasal cannula is dirty. 3. Cleaning detergent used to clean nasal cannulas may be absorbed into the plastic and can irritate the skin. 4. Nasal prongs are stiff causing nasal irritation and discomfort.	1. Contact your health care provider and/or doctor. 2. Wipe nasal cannula down with a damp cloth to remove oil and debris. If detergent is needed use a mild soap and rinse well. 3. Replace cannula. When cleaning cannula only use a damp cloth. Do not use strong detergents, disinfectants or oil based soaps. 4. Replace nasal cannula. Do not use a nasal cannula for more than 30 days.
Nasal prongs and tubing is stiff	1. Most nasal cannulas are made with a PVC material, which may harden with age and extended use. 2. Alcohol based cleaners may harden the PVC material	1. Replace your nasal cannula 2. Replace your nasal cannula

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4340-00		Power Supply: C8 (2-prong), 110-220V	for:	For Power Cord with C7 (2-Prong)	
		The 2-prong Power Supply is the standard power supply used with NanoVi and NanoVi Wellness devices			
4350-00 4350-10 4350-20 4350-30	Power Cord C7: to 2-prong 6.5 ft. (2.0m)  Type A Type C Type G Type I				
		Type A (US, CA, Mexico, Japan)	Type C (EU, Asia, Israel, South America)	Type G (UK, Malaysia, Singapore)	Type I (AUS, New Zealand)
4350-60	Power Cord C7: to 2-prong 6.5 ft. (2.0m)  Type A-CH		Type A-CH (China)		

## TR60M SERIES 60W MEDICAL SWITCHING ADAPTER

### Features

- \* Universal Input Range 90~264VAC
- \* Meets EN60601-1 and EN55011 Class B
- \* Continuous Short Circuit Protection
- \* Over Voltage Protection
- \* Meet CEC Level IV  
(Output Cable Length  $\leq$  1800mm)  
(TR60M Series meets CEC IV except TR60M05 is Non-CEC Compliant)  
(TR60M12: Output Cable Length  $\leq$  1220mm 16AWG)
- \* Efficiency & Standby Power Meet Level V (Option)  
(Output Cable Length  $\leq$  1800mm)  
(TR60M12 : Output Cable Length  $\leq$  720mm 16AWG)  
(TR60M15 : Output Cable Length  $\leq$  1220mm 16AWG)  
(TR60M18, TR60M19 : Output Cable Length  $\leq$  1500mm 18AWG)
- \* Meets 2MOPP

	36V	1.66A	360mVp-p	$\pm 2\%$	$\pm 1\%$	$\pm 2\%$	87%
TR60M48	48V	1.25A	480mVp-p	$\pm 2\%$	$\pm 1\%$	$\pm 2\%$	87%

**Specifications**

**INPUT SPECIFICATIONS:**

Voltage ..... 90~264Vac  
 Frequency ..... 47 to 63Hz  
 Inrush Current ..... Cold Start @25°C 80A max. @240Vac  
 Conducted EMI ..... CISPR/FCC Class B  
 Leakage Current ..... 0.1mA max.

**OUTPUT SPECIFICATIONS:**

Holdup Time ..... 8ms typ. @115Vac  
 Short Circuit Protection ..... Continuous  
 Over Voltage Protection ..... Yes  
 Temperature Coefficient ..... ±0.05%/°C


**GENERAL SPECIFICATIONS:**

Isolation ..... Input to output =5,656VDC  
 Operating Temperature ..... 0 ~ 60°C (see derating curve)  
 Storage Temperature ..... -20 ~ 85°C  
 Humidity ..... 93% RH max. Non condensing  
 Cooling ..... Natural Convection  
 Switching Frequency ..... 100KHz Typical  
 MTBF ..... MIL-HDBK-217F, GB, at 25°C/115VAC ..... 200K hrs min.  
 Altitude ..... 3000m  
 Dimensions ..... 5.197x2.283x1.201inches (132.00x58.00x30.50mm)  
 Weight ..... 345g(0.76 Pounds)

**SAFETY AND EMC:**

Emission and Immunity ..... EN55011, EN60601-1-2, EN61000-3-2  
 EN61000-3-3  
 Safety ..... IEC60601-1, EN60601-1, UL ANSI/AAMI ES60601-1:2005

**TR60M Series Derating Curve**


**NOTE:**

1. Voltage accuracy at 60% full load.
2. Add a 0.1uF ceramic capacitor and a 10uF E.L. capacitor to output for Ripple & Noise measurement @20MHz BW.
3. Line regulation measured from 100Vac to 240Vac, full load.
4. Load regulation measured from 60% to 100% full load and from 60% to 20% full load (60% +/- 40% full load).
5. Typical efficiency at 230VAC and full load at 25°C.
6. "Various TR Series adapters are PSE certified. PSE certification alone is not sufficient for importation into Japan. A valid PSE mark must contain the name of the importer as shown in the example below. If PSE mark is required, the name of the registered importer must be supplied to Cincon on order placement. Product labels will not contain PSE mark if importer name is not supplied. Consult factory or local representative for details".

**Mechanical Specification**

All Dimensions are in inches(mm)  
 Tolerance: Inches:X.XXX±0.02  
 Millimeters:X.XX±0.5

Typical at 25°C, nominal line and 75% load, unless otherwise Specified

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© Copyright 2025 Eng3 Corporation. All rights reserved. M126-rev11